

## **Cardiac Resynchronization Therapy: Evaluation and Management of Non- Responders**

**Paul A. Levine, M.D., FACC**  
Vice President, Medical Services  
St. Jude Medical CRMD  
Clinical Professor of Medicine  
Loma Linda University School of Medicine

Title slide for presentation on Pacemaker Diagnostics, a comparison of the various diagnostic features available from multiple different manufacturers. This series will NOT present absolutely every possible iteration and feature of every device as this is simply not feasible but it will provide an overview of the major features and capabilities of the larger manufacturers. Where possible, I have attempted to use recent printouts from the current devices but this will quickly become outdated.

These slides were all obtained from patients whom I have cared for. The commentary reflects the author's understanding of the various capabilities of these systems. Any inaccuracies are mine alone. These printouts are provided as a service to individuals with an interest in this field simply to increase an awareness of the various diagnostic capabilities that presently exist. Some of the commentary reflect my personal experience with the various diagnostic feature and counters. It is rare that a single diagnostic is adequate for the patient. Each has its own strengths and weaknesses - each tend to compliment the others and are best used in conjunction with one another to facilitate an understanding of the behavior of the pacing system in each patient.

Slide series 1998\98040706.ppt

Slide # PAL\980101

- Symptomatic congestive heart failure, NYHA Functional Class III-IV
- Maximal medical therapy
- End-diastolic LV dimension > 55 mm
- Intraventricular conduction defect (usually LBBB) with QRS > 130 ms

**In all the major studies and in general clinical practice 20-30% of patients do not respond to CRT**

## What is Cardiac Resynchronization?

- **Optimize the Left Atrial-Left Ventricular synchrony**
  - to maximize ventricular filling and eliminate diastolic regurgitation
- **Improve coordination of mechanical contraction of the Left Ventricle**
  - Stimulating LV septum (from RV) and LV posterior-lateral wall either simultaneously or sequentially
    - We will NEVER totally normalize LV contraction
  - CRT is NOT to coordinate RV and LV contraction
    - Mis-interpretation of “biventricular pacing”


- **Lack of mechanical dyssynchrony despite LBBB**
- **Poor LV lead location**
- **Large area of scar from prior M.I.**
  - Myocardium in the area of the LV lead must be able to contract
- **Pacemaker inhibited a significant portion of the time**
  - Atrial fibrillation with intact conduction
  - Automatic mode switch
  - Sinus rate > Maximum Tracking Rate
- **Inappropriate AV and VV delays**
- **LV lead dislodgment or loss of capture**

## Inappropriate Patient Selection

- While the ECG (LBBB) is a good screening tool, not all patients have mechanical dyssynchrony
- Large scar due to prior MI, remainder of myocardium may be hypercontractile
- Predominant diastolic rather than systolic dysfunction

**No options available for these patients; potential placebo effect (30% of control patients had beneficial effect)!**

- Lead placed in CS and advanced to anterior cardiac vein
  - This is biventricular but not CRT



## LV Lead problems

- **The majority of patients who receive a CRT system have an intact rhythm and AV nodal conduction**
  - This is likely to change in the future
- **Failure to deliver CRT will result in a return of the intrinsic conduction pattern**
  - While the etiology of a malfunction will be the same as other pacing system malfunctions, the symptomatic presentation will differ as heart rate will continue to be supported by the other ventricular lead
    - Some clinicians (Prof. J.J. Blanc) have suggested that LV pacing alone works as well as BiV
    - This is being evaluated by SJM in the B-Left Study in Europe

## Pacing System Malfunction

- **Changes in capture threshold**
  - Loss of capture
- **Alterations in sensing**
  - Undersensing
  - Oversensing
- **Lead dislodgment**
- **Timing considerations**
  - Paced and sensed AV delays
    - Optimize
  - Refractory periods
    - As short as possible
  - Maximum tracking rate
    - As high as possible

## Signs & Symptoms

### **Standard Pacing System**

- Syncope, Presyncope
- Fatigue, Low C.O.
- Pre-implant Symptoms
- Loss of Capture - asystole
- Pauses in paced rhythm

### **Biventricular Pacing System**

- Exacerbation of CHF
- Intact rhythm on ECG
- Classic paced complexes
- Return of intrinsic PR interval

**Pacing system malfunction in a CRT system may be subtle and difficult to identify from the ECG**


- **12 lead ECGs**
  - Native rhythm and conduction
  - RV paced rhythm
  - LV paced rhythm
  - Biventricular paced rhythm
- **Intracardiac electrograms**
- **Telemetered event markers**

## Loss of Capture



### **Standard pacing system**


- Pacing stimulus fails to elicit a cardiac depolarization



### **CRT System**


- Wider or subtle change in paced QRS morphology

## Loss of Capture: Utility of Event Markers and EGM




With current ICD-D systems, will only see this marker pattern with loss of RV capture since sensing has been disabled in the LV channel

## Ventricular Capture Threshold



**True biventricular AutoCapture is currently under-development.**

## Ventricular Capture Threshold



**With LV capture, there is conduction to the RV with "R" sense. Loss of LV capture results in the intrinsic rhythm.**

## Etiology of Loss of LV Capture

- **Lead dislodgment**
- **Massive increase in capture threshold**
- **Poor position requiring high outputs**
- **Mechanical problem with LV lead**
- **Pharmacologic / metabolic factors**


## Management of Loss of Capture

- **Increase output**
- **Correct metabolic abnormalities**
- **Withdraw recently introduced pharmacologic agents**
- **Program to unipolar output and sensing configuration**
  - Possible with CRT-P but not CRT-D
- **Operative intervention to reposition or replace lead**


- **Limiting the MTR is NOT the appropriate way to limit the highest heart rate**
  - IF the sinus rate can exceed MTR and there is intact AV nodal conduction, when the patient has a physiologic need for a high rate, there will be loss of CRT
  - In presence of AV block, sinus > MTR will cause a loss of appropriate AV synchrony even though Ventricular CRT is still present

## Timing: Maximum Tracking Rate

### Complete Heart Block




Sinus rate > MTR




## Sinus > Maximum Tracking Rate

### Congestive Heart Failure - CRT Pacing



Sinus rate > MTR; PV < PR interval



## Maximum Tracking Rate

Summary (Since Last Cleared)

Last Cleared: Aug 16, 2007 6:49 AM

Time Sampled	153d 6h 7m 32s
Percent of Counts Paced in Atrium	15%
Percent of Counts Paced in Ventricle	>99%
Total Time at Max Track Rate	1d 4h 33m 21s
AMS Episodes	1
PMT Detections (Passive)	0
RRP Intervals < 60 ms	0

P  
A  
C  
E  
D  
  
R  
A  
T  
E


MTR<sub>2</sub>

MTR

Base  
Rate

ATRIAL RATE

## Stored EGMs




**Native atrial rate > MTR**

**P wave in PVARP - hence not tracked but is conducted**

**In a CRT pacing system, set MTR >>  
maximum anticipated sinus rate**

## Timing: Assessment of AV Delay



Significant latency between atrial output and atrial depolarization

Atrial depolarization fuses with paced QRS

### Observations with respect to AV delay

First examine the complexes associated with loss of ventricular capture. The atrial output induces an atrial depolarization. However, there is significant latency between the atrial stimulus and the atrial evoked response. It measures almost 80 ms. At a programmed AV delay of 150 ms, that means that the P wave effectively starts 70 milliseconds in front of the paced QRS.


Now turn to the left side of the rhythm strip. The up-slope of the QRS complex that would normally be labeled an initial R wave is identified by a red arrow. This is really the first part of the P wave. The blue line aligning with the ventricular output marker identifies the start of the QRS which occurs after this "R" wave.

There is normally a delay between electrical activation of the atrium and mechanical contraction. In addition, with pacing from the Right Atrium, there is a further delay associated with the need to conduct the atrial depolarization from the RA to the LA. In all likelihood, at very short AV delays, the left atrial contraction will coincide with the LV contraction precluding any benefit from atrial transport. Although patients with a failing heart function on a flattened Starling or Ventricular Function curve, they commonly need every little bit of help that we can provide. In this case, it is to maintain optimal atrial transport to maximize ventricular filling.

## Programming AV / PV Delay Multisite Fusion


### AV Delay

100 ms



Fusion due to  
CRT Pacing

160 ms



Fusion due to  
CRT Pacing

210 ms




Fusion due to  
CRT Pacing  
and Native  
Conduction

Native



Intrinsic  
Conduction

## Programming AV / PV Delay using Echo-Doppler at PV 200 ms



### LV Inflow tract

VTI 0.27 m  
V 1.26 m/s



### LV Outflow tract

VTI 0.250 m

- Needs to be scheduled
- Time consuming
- 10-15% error between sequential measurements and between interpreters
- There is no consensus as to the best technique or measurement
  - Standard M mode (transthoracic or transesophageal)
  - Tissue Doppler Imaging (TDI)
  - Echo strain

**Although Echo-Doppler is the current gold-standard, it is not convenient and it cannot be done at each follow-up visit to adjust AV and VV intervals with reverse remodeling.**

- **QuickOpt is an algorithm incorporated in the SJM programmer to adjust paced and sensed AV delays and VV delays based on EGM measurements**
- **Contraction cannot occur until the muscle has been depolarized**
- **Measurement of the time to depolarization in each chamber from atrial depolarization and between chambers will provide a guide towards programming paced and sensed AV delays**

**In an analysis of 11 separately published studies, 440 out of 550 patients (80%) showed statistically significant improvement from sequential biventricular pacing over simultaneous pacing.<sup>1-11</sup>**

<sup>1</sup>Chan, et al. "Tissue Doppler Guided Optimization of A-V and V-V Delay of Biventricular Pacemaker Improves Response to Cardiac Resynchronization Therapy in Heart Failure Patients" *J of Cardiac Failure* 2004; 10, 4 (suppl.): S72 (abstract 199).

<sup>2</sup>Bordachar, et al. "Echocardiographic Parameters of Ventricular Dyssynchrony Validation in Patients with Heart Failure Using Sequential Biventricular Pacing" *JACC* 2004 Dec 7; 44 (11): 2157-2165.

<sup>3</sup>Vanderheyden, et al. "Tailored echocardiographic interventricular delay programming further optimizes left ventricular performance after cardiac resynchronization therapy" *Heart Rhythm*, Volume 2, No. 10, Oct 2005 1066-1072.

<sup>4</sup>Van Gelder, et al. "Effect of Optimizing the VV Interval on Left Ventricular Contractility in Cardiac Resynchronization Therapy" *American Journal of Cardiology* 2004; 93, 1500-1503.

<sup>5</sup>Sogaard, et al. "Sequential Versus Simultaneous Biventricular Resynchronization for Severe Heart Failure: Evaluation by Tissue Doppler Imaging" *Circulation* 106: 2078-2084 (2002).

<sup>6</sup>Rosanio, et al. "Non-Simultaneous Pacing of the Right and Left Ventricles for Heart Failure: Is It Worth It?" *AHA Abstract: 1618 (2003) AHA 76th Scientific Sessions*, Orlando, Nov. 9-12, 2003.

<sup>7</sup>Cochlain, et al. "The Effect of Variation in the Interval Between Right and Left Ventricular Activation on Paced QRS Duration" *PACE* 2001; 24: 1750-1752.

<sup>8</sup>Perego, et al. "Simultaneous vs. Sequential Biventricular Pacing in Dilated Cardiomyopathy: An Acute Hemodynamic Study" *The European Journal of Heart Failure* 2003; 5: 305-313.

<sup>9</sup>Bracke, et al. "Importance of Interventricular Delay to Optimize Cardiac Resynchronization Therapy" *JACC* 41: (2003) *ACC 52nd Annual Scientific Sessions March 30-April 2nd, 2003*, Chicago.

<sup>10</sup>Mortenson, et al. "Sequential Biventricular Pacing: Evaluation of Safety and Efficacy" *PACE* 2004; 27: 339-345.

<sup>11</sup>Leon, et al. "Effect of Cardiac Resynchronization Therapy with Sequential Biventricular Pacing on Doppler-Derived Left Ventricular Stroke Volume, Functional Status and Exercise Capacity in Patients with Ventricular Dysfunction and Conduction Delay" *PACE* 25: 141 (2002) *NASPE 23rd Annual Scientific Sessions, May 8-11, 2002*, San Diego.

- **Three studies comprising 350 patients studied patient quality of life (QoL) and all showed improvement on QoL scores<sup>1-3</sup>**
- **In addition, two studies specifically address non-responder rates and both showed a reduction in non-responder rates<sup>1,4</sup>**

<sup>1</sup>Chan, et al. "Tissue Doppler Guided Optimization of A-V and V-V Delay of Biventricular Pacemaker Improves Response to Cardiac Resynchronization Therapy in Heart Failure Patients" *J of Cardiac Failure* 2004; 10, 4 (suppl.): S72 (abstract 199).

<sup>2</sup>Bordachar, et al. "Echocardiographic Parameters of Ventricular Dyssynchrony Validation in Patients with Heart Failure Using Sequential Biventricular Pacing" *JACC* 2004 Dec 7; 44 (11): 2157-2165.

<sup>3</sup>Leon, et al. "Effect of Cardiac Resynchronization Therapy with Sequential Biventricular Pacing on Doppler-Derived Left Ventricular Stroke Volume, Functional Status and Exercise Capacity in Patients with Ventricular Dysfunction and Conduction Delay" *PACE* 25: 141 (2002) *NASPE 23rd Annual Scientific Sessions, May 8-11, 2002, San Diego*.

<sup>4</sup>Vanderheyden, et al. "Tailored echocardiographic interventricular delay programming further optimizes left ventricular performance after cardiac resynchronization therapy" *Heart Rhythm*, Volume 2, No. 10, Oct 2005 1066-1072.

- Seven studies used NYHA class improvement as a surrogate for patient response to optimized therapy
- All seven studies resulted in improvements in NYHA functional class status<sup>1-7</sup>

<sup>1</sup>Chan, et al. "Tissue Doppler Guided Optimization of A-V and V-V Delay of Biventricular Pacemaker Improves Response to Cardiac Resynchronization Therapy in Heart Failure Patients" *J of Cardiac Failure* 2004; 10, 4 (suppl.): S72 (abstract 199).

<sup>2</sup>Bordachar, et al. "Echocardiographic Parameters of Ventricular Dyssynchrony Validation in Patients with Heart Failure Using Sequential Biventricular Pacing" *JACC* 2004 Dec 7; 44 (11): 2157-2165.

<sup>3</sup>Vanderheyden, et al. "Tailored echocardiographic interventricular delay programming further optimizes left ventricular performance after cardiac resynchronization therapy" *Heart Rhythm*, Volume 2, No. 10, Oct 2005 1066-1072.

<sup>4</sup>Sogaard, et al. "Sequential Versus Simultaneous Biventricular Resynchronization for Severe Heart Failure: Evaluation by Tissue Doppler Imaging" *Circulation* 106: 2078-2084 (2002).

<sup>5</sup>Rosano, et al. "Non-Simultaneous Pacing of the Right and Left Ventricles for Heart Failure: Is It Worth It?" *AHA Abstract: 1618* (2003) *AHA 76th Scientific Sessions, Orlando, Nov. 9-12, 2003*.

<sup>6</sup>Mortenson, et al. "Sequential Biventricular Pacing: Evaluation of Safety and Efficacy" *PACE* 2004; 27: 339-345.


<sup>7</sup>Leon, et al. "Effect of Cardiac Resynchronization Therapy with Sequential Biventricular Pacing on Doppler-Derived Left Ventricular Stroke Volume, Functional Status and Exercise Capacity in Patients with Ventricular Dysfunction and Conduction Delay" *PACE* 25: 141 (2002) *NASPE 23rd Annual Scientific Sessions, May 8-11, 2002, San Diego*.

## Timing Cycle Optimization

### Predicated upon

- Intact LV capture
- Proper device function

Timing cycle optimization is the first line of defense for non-responders and also improves the outcomes of patients that do respond to CRT therapy, particularly as reverse remodeling occurs.



Aranda, Juan M. ISHNE Heart Failure World-Wide Internet Symposium, April 2006.


ISHNE Heart Failure World-Wide Internet Symposium and Noninvasive Electrocardiology

Assuming that LV capture is maintained and the device is functioning properly, QuickOpt is your opportunity to convert non-responders and to improve the outcomes of patients who do respond to CRT therapy

## Timing Cycle Optimization

### Optimal Delays Change Acutely...and Often<sup>1</sup>

- 63 pts, EF < 35%
- NYHA ≥ II, QRS > 150 ms
- LV lead in lateral or postero-lateral vein
- Results:
  - Only 3 pts unchanged
  - 18 pts needed adjustments at each FU
  - VV 73 times in 27 pts
  - AV 43 times in 21 pts



<sup>1</sup>O'Donnell, et al. "Long-Term Variations in Optimal Programming of Cardiac Resynchronization Therapy Devices" PACE: Vol 28 Supp S24-S26 (Jan 2005).

ISHNE International Society for Holter and Noninvasive Electrocardiology

31

In this study by O'Donnell, appearing in PACE in 2005, you can see that the question of whether timing cycles change over time is answered.

In some patients, the timing cycles began to change in a little as 24 hours and as often as every two to three weeks. As you can see, the VV delays appear to get shorter while the AV delays tend to get longer.

In this group of 63 patients, only three (less than 5%) of the patients were unchanged. 18 patients needed adjustment at every follow up with VV delays being changed 73 times in 27 patients and AV delays being changed 43 times in 21 patients over a 9 month period.

## QuickOpt™ Timing Cycle Optimization

- **QuickOpt™ optimization is clinically-proven to correlate with echo based methods:**
  - Prospective 11-patient pilot study presented at Cardiosim 2004<sup>1</sup>
  - Retrospective 61-patient study conducted in the U.S.<sup>2</sup>
  - IDE clinical trial with 115 patients conducted in the U.S.<sup>3</sup>

<sup>1</sup> Meine, et al. "IEGM Based Method for Estimating Optimal VV Delay in Cardiac Resynchronization Therapy." *Europace Supplements*, Vol. 6, June 2004 (#149/2).

<sup>2</sup> Meine, et al. "An Intracardiac EGM Method for VV Optimization During Cardiac Resynchronization Therapy" *Heart Rhythm Journal* 3 (5) May 2006 [abstract AB30-5].

<sup>3</sup> Porterfield, et al. "Device based intracardiac delay optimization vs. echo in ICD patients (Acute IEGM AV/PV and VV Study)" *Europace Vol 8 Supp 1* July 2006 [abstract #6178].

### Study Method<sup>1-3</sup>


- **Echo-derived Optimization**
  - Aortic VTI measured at several AV and VV delay settings
  - Optimal delay determined by the maximum aortic VTI (Max AVTI)
- **SJM IEGM-derived VV Optimization**
  - SJM IEGM formula used to determine optimal AV and VV delay
  - Aortic VTI measured at IEGM-derived delay
- **Aortic VTI at IEGM-derived VV delay compared to Max AVTI from echo optimization**
- **Correlation measures how closely IEGM VTI matches echo VTI**
- **Aortic VTI measures forward flow!**

<sup>1</sup> Meine, et al. "IEGM Based Method for Estimating Optimal VV Delay in Cardiac Resynchronization Therapy." *Europace Supplements*, Vol. 6, June 2004 (#149/2).

<sup>2</sup> Meine, et al. "An Intracardiac EGM Method for VV Optimization During Cardiac Resynchronization Therapy" *Heart Rhythm Journal* 3 (5) May 2006 [abstract AB30-5].


<sup>3</sup> Porterfield, et al. "Device based intracardiac delay optimization vs. echo in ICD patients (Acute IEGM AV/PV and VV Study)" *Europace* Vol 8 Supp 1 July 2006 [abstract #6178].

### A Pilot Prospective Study - Cardiostim 2004<sup>1</sup>



<sup>1</sup>Meine, et al. "IEGM Based Method for Estimating Optimal VV Delay in Cardiac Resynchronization Therapy." Europace Supplements, Vol. 6, June 2004 (#149/2).

## Results of Retrospective IEGM VV Optimization Study<sup>1</sup>




<sup>1</sup>Meine, et al. "An Intracardiac EGM Method for VV Optimization During Cardiac Resynchronization Therapy" Heart Rhythm Journal 3 (5) May 2006 [abstract AB30-5].

## QuickOpt™ Optimization


### ACUTE IEGM AV and PV (Sensed AV) Study<sup>1</sup>

- PV Delay Group (n=56)
  - CCC = 96.1%
- AV Delay Group (n=56)
  - CCC = 97.5%

Plot of Max Echo AVTI for the PV Delay on 52 Analyzable Patients with the Identity Line Superimposed



Plot of Max Echo AVTI for the PV Delay on 52 Analyzable Patients with the Identity Line Superimposed




<sup>1</sup>Porterfield, et al. "Device based intracardiac delay optimization vs. echo in ICD patients (Acute IEGM AV/PV and VV Study)" Europace Vol 8 Supp 1 July 2006 [abstract #6178].

## ACUTE IEGM-CRT VV Study<sup>1</sup>

- VV Delay Group (n=54)
  - CCC = 96.6

Plot of Max Echo AVTI for the PV Delay on 52 Analyzable Patients with the Identity Line Superimposed



<sup>1</sup>Porterfield, et al. "Device based intracardiac delay optimization vs. echo in ICD patients (Acute IEGM AV/PV and VV Study)" Europace Vol 8 Supp 1 July 2006 [abstract #6178].

## QuickOpt™ Timing Cycle Optimization

- **QuickOpt™ optimization is clinically-proven to correlate with echo based methods:**

- **Prospective 11-patient pilot study<sup>1</sup>**
  - 99.2% correlation
- **Retrospective 61-patient study<sup>2</sup>**
  - 97.6% correlation
- **Prospective multi-center IDE clinical trial<sup>3</sup>**
  - 96.1% correlation for Sensed AV (PV) delay<sup>3</sup>
  - 97.5% correlation for Paced AV delay<sup>3</sup>
  - 96.6% correlation for VV Delay<sup>3</sup>

**Special Caution: Leads in RV Apex and Rt. Atrial Appendage, other lead locations have not yet been studied!**

<sup>1</sup> Meine, et al. "IEGM Based Method for Estimating Optimal VV Delay in Cardiac Resynchronization Therapy." *Europace Supplements*, Vol. 6, June 2004 (#149/2).

<sup>2</sup> Meine, et al. "An Intracardiac EGM Method for VV Optimization During Cardiac Resynchronization Therapy" *Heart Rhythm Journal* 3 (5) May 2006 [abstract AB30-5].


<sup>3</sup> Porterfield, et al. "Device based intracardiac delay optimization vs. echo in patients (Acute IEGM AV/PV and VV Study)" *Europace Vol 8 Supp 2006* [abstract #6178].




**What if non-invasive adjustments do not  
improve response rate?**

- Detailed echo study (TDI, strain, speckle) looking for area of latest contraction
- Placement of a second lead in the area of latest contraction
  - Transvenous if possible
  - Epicardial

## Addition of a second LV lead via the CS



**Baseline**



**Following revision**

Courtesy – George Eisinger, FCE – SJM (2007)

## **Summary: CRT Non-responders**

- **If patient presents with exacerbation of CHF - pacing system malfunction is now part of differential diagnosis**
- **Baseline 12 lead ECGs provide valuable templates for initial screening**
- **Change in QRS morphology or inhibition of pacing system suggests a problem necessitating further evaluation and may allow a return of CHF symptoms**

### Summary (continued): CRT non-responders


- For patients who do not respond to CRT or who, after a period of response, deteriorate
- Keep refractory periods as short as possible and MTR as high as possible
- Optimize AV and VV intervals
  - Echo-Doppler
  - QuickOpt (can be done at each visit)



International Society for Holter  
and Noninvasive Electrocardiology

## Supplemental Slides on QuickOpt™

## Electrical Conduction Pathway



## Key Principles

- Mitral valve closure can be estimated by measuring the interatrial conduction time (P-wave duration)<sup>1</sup>
- Onset of isovolumic contraction can be measured using the peak of the R wave<sup>2</sup>
- Interventricular conduction delays can be measured by evaluating simultaneous RV & LV IEGMs and measuring the time between the peaks of the R waves<sup>3</sup>

<sup>1</sup>Worley, et.al "Optimization of cardiac resynchronization: left atrial electrograms measured at implant eliminates the need for echo and identifies patients where AV optimization is not possible" *Journal of Cardiac Failure Aug. 2004 Vol. 10, Issue 4, Pg S62.*

<sup>2</sup>Koglek, et al. "Eine einfache methode zur bestimmung des AV-intervalls bei zweikammerschrittmachern," *Herzschrittmacherther Elektrophysiol* 2000; 11:244-253.

<sup>3</sup>Meine, et al. "IEGM Based Method For Estimating Optimal VV Delay in Cardiac Resynchronization Therapy" *Europace Supplements*, Vol. 6, June 2004 (#149/2).

- **Characterize interatrial conduction patterns in order to:**
  - Maximize preload
  - Allow proper timing of mitral valve closure
- **Characterize intrinsic and paced interventricular conduction patterns so that:**
  - Pacing stimuli and the resultant LV & RV conduction (paced wave fronts) meet midway between the two stimulation sites

## Electronic Optimization of Sensed AV Delays

The IEGM estimates inter-atrial conduction delay by measuring the width of the P wave. The QuickOpt™ algorithm utilizes this measurement to calculate the optimal sensed and paced AV delays, with the goal maximizing preload and allowing for proper timing of mitral valve closure.

The optimal sensed AV delay is the sum of  $A_s$  (sensed P wave) and  $\Delta$  (an added interval that allows for ventricular filling). The interval added is from 30 to 60 ms, depending on the measured P Wave. A short P wave would have a longer interval added, while a longer P wave would have a shorter interval added.


$$SAV_{opt} = A_s + \Delta$$

$$30 \leq \Delta \leq 60$$


SAV is Sensed AV Delay

$A_s$  is sensed P wave

$\Delta$  is added interval




## Electronic Optimization of Sensed AV Delays




### Short P Wave Example

$$A_S + \Delta = SAV_{opt}$$
$$50 + 60 = 110$$



### Long P Wave Example

$$A_S + \Delta = SAV_{opt}$$
$$110 + 30 = 140$$



$$SAV_{opt} = A_S + \Delta$$
$$30 \leq \Delta \leq 60$$

## Electronic Optimization of Paced AV Delays

The optimal paced AV delay is the sum of  $A_S$  (sensed P wave) and  $\Delta$  (an added interval that allows for ventricular filling). The interval added is from 80 to 110 ms, depending on the measured P Wave. A short P wave would have a longer interval added, while a longer P wave would have a slightly shorter interval added.

$$PAV_{opt} = A_S + \Delta$$


$$80 \leq \Delta \leq 100$$

PAV is Paced AV Delay  
 $A_S$  is sensed P wave  
 $\Delta$  is added interval

### Short P Wave Example

$$A_S + \Delta = PAV_{opt}$$


$$50 + 110 = 160$$



### Long P Wave Example

$$A_S + \Delta = PAV_{opt}$$

$$110 + 80 = 190$$



$$PAV_{opt} = A_S + \Delta$$

$$80 \leq \Delta \leq 110$$

## V-V optimization


Paced and sensed tests are performed to characterize the conduction properties of the ventricles. The goal is to time the right and left ventricular activation so that the paced wavefronts meet near the ventricular septum. The V Sense test measures the intrinsic interventricular delay ( $\Delta$ ). The RV Pace and LV Pace tests are used to measure the difference in interventricular conduction delays ( $\epsilon$ ). All of these measurements are used to optimize the VV Delay (Interventricular Pace Delay).

$$VV_{\text{opt}} = 0.5 \times (\Delta + \epsilon)$$


$\Delta$  is related to the intrinsic depolarization

$\epsilon$  is a correction term depending on wave front velocity


## V-V optimization: Intrinsic depolarization term ( $\Delta$ )



## V-V optimization: wave front velocities ( $\varepsilon$ )



## V-V optimization: wave front velocities ( $\varepsilon$ )



$$VV_{opt} = 0.5 \times (\Delta + \varepsilon)$$

$\Delta$  is related to the intrinsic depolarization

$\varepsilon$  is a correction term depending on wave front velocity

**Example:**

$\Delta=40\text{ms}$  (RV sensed before LV)

IVCD\_LR=60ms      IVCD\_RL=50ms

In this case, it takes longer to depolarize from Left to Right


$\varepsilon = \text{IVCD\_LR} - \text{IVCD\_RL} = 60 - 50 = 10 \text{ ms}$

$VV_{opt} = 0.5 \times (\Delta + \varepsilon) = 0.5 \times (40 + 10) = 50/2 = 25$

We get: LV first, V-V=25ms

## How QuickOpt™ Optimization Works

QuickOpt™ Optimization electrically characterizes the conduction properties of the heart and uses an innovative formula to calculate the optimal AV and VV delays.



## How QuickOpt™ Optimization Works

### **AV Optimization:**

1. QuickOpt™ optimization measures the total P-wave duration of eight IEGM events for the A-Sense test
2. Measured P-wave durations are averaged
3. The QuickOpt™ algorithm uses an innovative formula to calculate the optimal AV and PV delays
  - A variable delta is added to the measured P-wave duration
  - Delta added is dependent on the measured events

## How QuickOpt™ Optimization Works

### **VV Optimization:**

1. QuickOpt™ optimization measures eight IEGM events for each of the V Sense, RV Pace and LV Pace tests.
  - V Sense – measures intrinsic interventricular delay
  - RV Pace – measures conduction speed from right to left
  - LV Pace – measures conduction speed from left to right
2. Measurements from each test are averaged
3. The QuickOpt™ algorithm uses an innovative formula to calculate the optimal VV delay
  - Difference in conduction speeds between RV Pace and LV Pace tests used as correction term for the intrinsic delay measurement

## QuickOpt™ – Brady Parameters Screen

QuickOpt™ Optimization is found on the brady parameters screen

Brady Tachy Special Episode Settings 70 bpm 852 ms

Mode: DDD Sensor: Off

**Pacing Rates**  
Base Rate 60 ppm / 1000 ms  
Hysteresis Rate Off  
Rest Rate Off  
Max Track Rate 90 ppm / 667 ms  
2:1 Block Rate 110 bpm / 545 ms

**Delays**  
AV/PV Delay 140 ms / 110 ms  
Rate Responsive AV/PV Delay Off  
Negative AV/PV Hysteresis -20 ms (>= 100 ms)

**Pacing Output**  
A. Output 2.5 V, 0.5 ms  
V. Output RV: 2.5 V, 0.6 ms / LV: 2.5 V, 0.5 ms  
V. Sensing RV Only  
V. Pacing Biventricular  
Interventricular Pace Delay Simultaneous  
LV Pulse Configuration LV Bipolar

**Extended**  
Auto Mode Switch DOI  
Atrial Tachycardia Detection Rate 225 bpm / 267 ms  
PVC Options Off  
PMT Options Passive  
PMT Detection Rate 90 bpm / 667 ms  
Ventricular Noise Reversion Mode VOO  
Ventricular Safety Standby On

**Post-Shock Pacing**  
Post-Shock Mode Off  
Post-Shock Pause 2 sec

**Refractory**  
Ventricular Pace Refractory 400 ms  
PVARP 440 ms  
Rate Responsive PVARP/VPREF Off

Interrogate Temporary Pacing **QuickOpt™ Optimization** Program

Summary Parameters Diagnostics Episodes Tests Measured Data Main Menu End Session

## QuickOpt™ – Start Test Screen

Tests can be run automatically or manually  
Parameters required for testing are adjusted automatically

The screenshot displays the 'QuickOpt™ Timing Cycle Optimization' dialog box. At the top, there are tabs for 'Brady', 'Tachy', 'Special', and 'Episode Settings'. Below these, 'Mode' is set to 'DDD' and 'Sensor' is 'Off'. A status indicator in the top right corner shows '70 bpm' and '852 ms'. The dialog box contains the following text:

To optimize AV Delay, PV Delay, and Interventricular Pace Delay, five sets of measurements need to be obtained.

The following parameter values will be used temporarily to obtain proper rhythms:

Mode: DDD  
Base Rate: 90 bpm and 40 bpm  
AV/PV Delay: 300 ms and 100 ms

If these values do not succeed, you will be prompted to select other values. To begin, press "Start Test". The test can be cancelled at any time.

At the bottom of the dialog box, there are four buttons: 'Start Test' (highlighted with a red circle), 'View Test Results', 'Run Manually' (also highlighted with a red circle), and 'Cancel'. The background of the screen shows various settings for 'Pacemaker' and 'Sensors'.

## QuickOpt™ – Test in Progress

Automatic testing performs all relevant tests

QuickOpt™ Timing Cycle Optimization

Now performing measurements...

- ✓ A Sense: Mode DDD, Base Rate 40 bpm, AV/PV Delays 300 ms
- ✓ A Pace: Mode DDD, Base Rate 90 bpm, AV/PV Delays 300 ms
- V Sense: Mode DDD, Base Rate 60 bpm, AV/PV Delays 300 ms
- RV Pace: Mode DDD, Base Rate 60 bpm, AV/PV Delays 100 ms
- LV Pace: Mode DDD, Base Rate 60 bpm, AV/PV Delays 100 ms

40%

Skip Measurement Cancel

## QuickOpt™ – Results

Results are displayed in about a minute<sup>1</sup>

QuickOpt™ Timing Cycle Optimization: Results

QuickOpt™ optimization measurements were successful!

	Programmed:	Optimal:
<input checked="" type="checkbox"/> AV Delay:	150 ms	170 ms
<input checked="" type="checkbox"/> PV Delay:	100 ms	120 ms
<input checked="" type="checkbox"/> Interventricular Pace Delay:	Simultaneous	15 ms (LV First)

Press "Select Values" to select optimized values for programming. These can be adjusted before programming.

Summary  
Parameters  
Diagnostics  
Episodes  
Tests  
Measured Data  
Main Menu  
End Session

<sup>1</sup>AV Optimization in dual-chamber devices take approximately 30 seconds.  
AV and VV Optimization in CRT-D devices take approximately 90 seconds.

## QuickOpt™ – Details

Test details and electrograms can be viewed individually and tests can be run manually

QuickOpt™ Timing Cycle Optimization: Details

70 bpm  
855 ms

AV/PV Delay Optimization

Atrial Wave Duration Tests

A Sense: 61 ms (62, 62, 62, 62, 62, 55, 62 ms)		AV Delay: 150 ms	Optimize: 170 ms
A Pace: EGM stored		PV Delay: 100 ms	Optimize: 120 ms

Interventricular Pace Delay Optimization

RV-LV Interval Tests

V Sense: 22 ms (LV First) (23, 23, 23, 31, 23, 23, 8, 23 ms)		Interventricular Pace Delay: Simultaneous	Optimize: 16 ms (LV First)
RV Pace: 16 ms (16, 16, 16, 16, 16, 16, 16 ms)			
LV Pace: 4 ms (0, 16, 0, 0, 0, 23, 0, 0 ms)			

Print Report

Close

Summary  
Parameters  
Diagnostics  
Epicodes  
Tests  
Measured Data  
Main Menu  
End Session

## QuickOpt™ – Program

Once values are determined, programming is easy

Brady		Tachy		Special		Episode Settings	
Mode: <b>DDD</b>		Sensor: <b>Off</b>				70 bpm 859 ms	
<b>Pacing Rates</b> Base Rate: 60 ppm / 1000 ms Hysteresis Rate: Off Rest Rate: Off Max Track Rate: 120 ppm / 500 ms 2:1 Block Rate: 149 bpm / 403 ms <b>Delays</b> AV/PV Delay: 170 ms / 120 ms Rate Responsive AV/PV Delay: Off Negative AV/PV Hysteresis: Off		<b>Extended</b> Auto Mode Switch: DDI Atrial Tachycardia Detection Rate: 225 bpm / 267 ms PVC Options: Off PMT Options: A Pace on PMT PMT Detection Rate: 120 bpm / 500 ms Ventricular Noise Reversion Mode: Pacer Off Ventricular Safety Standby: On				<b>Summary</b> <b>Parameters</b> <b>Diagnostics</b> <b>Episodes</b> <b>Tests</b> <b>Measured Data</b>	
<b>Pacing Output</b> A. Output: 2.5 V, 0.5 ms V. Output: RV: 2.0 V, 0.5 ms / LV: 4.0 V, 1.0 ms V. Sensing: RV Only V. Pacing: Biventricular Interventricular Pace Delay: 15 ms (LV First) LV Pulse Configuration: LV Bipolar <b>Refractory</b> Ventricular Pace Refractory: 250 ms PVARP: 280 ms Rate Responsive PVARP/VREF: Off		<b>Post-Shock Pacing</b> Post-Shock Mode: VVI Post-Shock Base Rate: 75 ppm / 800 ms Post-Shock Pause / Duration: 1 sec / 30 sec Post-Shock V. Output: 5.0 V / 1.0 ms				<b>Main Menu</b> <b>End Session</b>	
<input type="checkbox"/> /U/L Interrogate		<input type="button" value="Clear Selected"/>		<input type="button" value="Temporary Pacing"/>		<input type="button" value="QuickOpt™ Optimization"/>	
		<input type="button" value="Program Preview"/>		<input checked="" type="button" value="Program"/>			

International Society for Holter and Noninvasive Electrocardiology

Additional Data  
Supporting the Effectiveness  
of Timing Cycle Optimization

## **Additional Data Supporting the Effectiveness of Timing Cycle Optimization**

- 45 patients classified as “non-responders” following implant**
- A-V & V-V delay optimized using Echo TDI**
- 85% of patients improved significantly after optimization as shown by QoL and EF**

Chan et al. “Tissue Doppler Guided Optimization of A-V and V-V Delay of Biventricular Pacemaker Improves Response to Cardiac Resynchronization Therapy in Heart Failure Patients” J of Cardiac Failure 2004; 10, 4 (suppl.): S72 (abstract 199).

## Additional Data Supporting the Effectiveness of Timing Cycle Optimization

- 41 HF patients with CRT device implantation
- Echo w/ TT & pulsed Doppler tissue imaging performed
  - Pre-CRT implantation
  - Post-CRT implantation w/ AV Optimization
  - 3 months post AV & VV optimization
- Results (After 3 mos. Optimized Therapy)
  - NYHA functional class improved ( $1.7 \pm 0.6$  vs.  $3.2 \pm 0.5$   $p < 0.01$ )
  - QoL scores Improved ( $24 \pm 14$  vs.  $45 \pm 17$   $p < 0.01$ )
  - Exercise Capacity Improved ( $389\text{m} \pm 65\text{m}$  vs.  $264\text{m} \pm 73\text{m}$   $p < 0.05$ )
  - LVEF Increased ( $34\% \pm 6\%$  vs.  $28\% \pm 6\%$   $p = 0.035$ )
  - 85% of patients were optimized at a V-V delay other than zero
- Individually optimized sequential CRT provided a significant hemodynamic improvement compared with simultaneous CRT

Bordachar et al. "Echocardiographic Parameters of Ventricular Dyssynchrony Validation in Patients with Heart Failure Using Sequential Biventricular Pacing" *JACC* 2004 Dec 7; 44 (11): 2157-2165.



## Additional Data Supporting the Effectiveness of Timing Cycle Optimization

	Simultaneous CRT	Optimized V-V	% Improvement	P value
VTI (mm)	122 ± 31	154 ± 42	26%	<.001
LVFT (ms)	404 ± 102	472 ± 110	17%	.001
LVFT/HR (ms/bpm)	6.0 ± 2.0	7.1 ± 2.2	18%	.001
IVD	35 ± 33	13 ± 25	63%	.013
VD	51 ± 34	34 ± 18	33%	.010
SUM dyssynchrony (ms)	86 ± 49	47 ± 31	45%	.002

Non-responder rate at 6 months = 10%

Vanderheyden, et al. "Tailored echocardiographic interventricular delay programming further optimizes left ventricular performance after cardiac resynchronization therapy" *Heart Rhythm*, Volume 2, No. 10, Oct 2005 1066-1072.



## Additional Data Supporting the Effectiveness of Timing Cycle Optimization

### 53 patients

#### 41 patients in sinus rhythm

26 patients w/ischemic cardiomyopathy (IC)

15 patients w/idiopathic dilated cardiomyopathy (IDC)

#### 12 patients in AF (23%)

#### Maximum dP/dt measured

Baseline (prior to CRT implant)

Simultaneous CRT pacing

V-V Optimized

	n	Baseline	Simultaneous	% Increase	Optimized V-V	% Increase
AF	12	941 ± 240	1,142 ± 207	21	1,180 ± 196	25
Sinus Rhythm	41	814 ± 178	952 ± 271	17	1,027 ± 297	26
IDC	15	754 ± 220	890 ± 245	18	955 ± 267	27
IC	26	846 ± 249	987 ± 232	17	1,069 ± 252	26
All Patients	53	842 ± 231	995 ± 247	18	1,061 ± 259	26

Van Gelder, et al. "Effect of Optimizing the VV Interval on Left Ventricular Contractility in Cardiac Resynchronization Therapy" *Am J Cardiol* 2004; 93, 1500-1503



## Additional Data Supporting the Effectiveness of Timing Cycle Optimization

In the van Gelder study, Simultaneous-CRT was compared to Sequential-CRT measuring LV dP/dt as the acute endpoint<sup>1</sup>

Rhythm	Disease	Baseline	Simul- taneous	Sequential	Relative Improve- ment	Simul- taneous	LV First	RV First
Sinus 15 pts	DCM	754±220	18%	27%	48%	4	10	1
Sinus 26 pts	ICM	846±249	17%	26%	58%	1	25	0
AF 12 pts		941±240	21%	25%	17%	2	9	1

87% of all patients were optimized at a V-V delay other than simultaneous.

<sup>1</sup>van Gelder, et al. Effect of optimizing the V-V interval on LV contractility in cardiac resynchronization therapy. Am. J Cardiology 2004, vol. 93: 1500-1503.

## Additional Data Supporting the Effectiveness of Timing Cycle Optimization

- **Both Studies:**


- Single center studies of consecutive patients with NYHA III/IV HF, QRS >130
- Trans-mitral flow optimized with AV delay post implant

- **Rosanio, et al.** AHA Abstract; 1618 (2003)

- N=22
- V-V Delay of 0 for first 2 months
- Echo based optimization of V-V delay (OPT) at 2 months

- **Sogaard, et al.** Circ. 2002; 106:2078-84

- N=20
- Optimal V-V Delay based on Echo TDI
- Acute data shown. After 3 mo. LVEF further improved to 39% (P<0.01)



## Additional Data Supporting the Effectiveness of Timing Cycle Optimization

**Bordachar, Rosanio and Sogaard studied SQ-CRT comparing  
EF, NYHA and 6 MWT<sup>1-3</sup>**

Author	# pts	Baseline			3 months			Optimal Sequence of Activation		
		% EF	NYHA	6 MWT	% EF	NYHA	6 MWT	LV	RV	Sim
Bordachar Improvement	41	28%	3.2	264 m	34% 21.4%	1.7 46.9%	369 m 39.8%	25	6	10
Rosanio Improvement	22	22.5 %	3.4	349 m	35% 55.5%	1.9 44.1%	441 m 26.4%	18	1	3
Sogaard Improvement	20	22.4 %	3.45	222 m	38% 30%	1.9 45%	401 m 80%	9	11	

**84% of all patients were optimized at a V-V delay other than simultaneous**



Orange Lacewing  
Feb 12, 2008  
Butterfly Farm, St. Thomas, USVI  
Native of Australia